

Conversión de potencia de LIOs del saco al sulcus



Power at Capsular Bag	Power at Ciliary Sulcus	Bag / Sulcus Difference	Subtract from Bag Power
+30.00 D	+28.55 D	-1.45 D	-1.50 D
+29.50 D	+28.09 D	-1.42 D	-1.50 D
+29.00 D	+27.61 D	-1.39 D	-1.50 D
+28.50 D	+27.14 D	-1.36 D	-1.50 D
+28.00 D	+26.67 D	-1.33 D	-1.00 D
+27.50 D	+26.20 D	-1.30 D	-1.00 D
+27.00 D	+25.73 D	-1.27 D	-1.00 D
+26.50 D	+25.26 D	-1.25 D	-1.00 D
+26.00 D	+24.79 D	-1.22 D	-1.00 D
+25.50 D	+24.31 D	-1.19 D	-1.00 D
+25.00 D	+23.84 D	-1.16 D	-1.00 D
+24.50 D	+23.36 D	-1.13 D	-1.00 D
+24.00 D	+22.89 D	-1.11 D	-1.00 D
+23.50 D	+22.42 D	-1.08 D	-1.00 D
+23.00 D	+21.94 D	-1.05 D	-1.00 D
+22.50 D	+21.47 D	-1.03 D	-1.00 D
+22.00 D	+21.00 D	-1.00 D	-1.00 D
+21.50 D	+20.53 D	-0.97 D	-1.00 D
+21.00 D	+20.05 D	-0.95 D	-1.00 D
+20.50 D	+19.58 D	-0.92 D	-1.00 D
+20.00 D	+19.11 D	-0.89 D	-1.00 D
+19.50 D	+18.63 D	-0.87 D	-1.00 D
+19.00 D	+18.16 D	-0.84 D	-1.00 D
+18.50 D	+17.69 D	-0.82 D	-1.00 D
+18.00 D	+17.21 D	-0.79 D	-1.00 D
+17.50 D	+16.73 D	-0.77 D	-1.00 D
+17.00 D	+16.26 D	-0.74 D	-0.50 D
+16.50 D	+15.78 D	-0.72 D	-0.50 D
+16.00 D	+15.31 D	-0.69 D	-0.50 D
+15.50 D	+14.83 D	-0.67 D	-0.50 D
+15.00 D	+14.35 D	-0.64 D	-0.50 D
+14.50 D	+13.88 D	-0.62 D	-0.50 D
+14.00 D	+13.40 D	-0.60 D	-0.50 D
+13.50 D	+12.93 D	-0.57 D	-0.50 D
+13.00 D	+12.45 D	-0.55 D	-0.50 D
+12.50 D	+11.97 D	-0.53 D	-0.50 D
+12.00 D	+11.49 D	-0.50 D	-0.50 D
+11.50 D	+11.02 D	-0.48 D	-0.50 D
+11.00 D	+10.54 D	-0.46 D	-0.50 D
+10.50 D	+10.07 D	-0.43 D	-0.50 D
+10.00 D	+9.58 D	-0.41 D	-0.50 D
+9.50 D	+9.11 D	-0.39 D	-0.50 D
+9.00 D	+8.63 D	-0.37 D	No Change
+8.50 D	+8.16 D	-0.35 D	No Change
+8.00 D	+7.68 D	-0.32 D	No Change
+7.50 D	+7.20 D	-0.30 D	No Change
+7.00 D	+6.72 D	-0.28 D	No Change
+6.50 D	+6.24 D	-0.26 D	No Change
+6.00 D	+5.76 D	-0.24 D	No Change
+5.50 D	+5.28 D	-0.22 D	No Change
+5.00 D	+4.81 D	-0.20 D	No Change